


Step 1 - UltraCheck™ 1 Calibration

Perform one UltraCheck 1 calibration per day or per each set of samples analyzed.




NOTE: If $RLU_{ATP1} \leq 5,000$ using a PhotonMaster or Lumitester C-110, rehydrate a new bottle of Luminase for maximum sensitivity.

Step 2 – Total ATP (tATP™) Analysis (1 per sample)

Included in QG21I™ and QG21It™ test kits.


2.1 – EXTRACTION

Add sample to extract ATP.



2.2 – DILUTION

Dilute out interferences.



2.3 – ASSAY

Measure ATP concentration.



NOTE: If $RLU_{dATP} \leq 10$ using a PhotonMaster or Lumitester C-110, you are below the low- detection limit.

NOTE: If $RLU_{tATP} \leq 50$ using a PhotonMaster or Lumitester C-110, consider accounting for background (RLU_{bg}). See Test Kit Instructions for guidance.

Total ATP (tATP) Calculation:


$$tATP (pg \text{ ATP/mL}) = \frac{RLU_{tATP}}{RLU_{ATP1}} \times 10,000 (pg \text{ ATP / mL})$$

Step 3 – Dissolved ATP (dATP™) Analysis (1 per sample)

Included in QG21I™ test kit only.


3.1 – DILUTION

Add sample to recover ATP.



3.2 – ASSAY

Measure ATP concentration.



NOTE: If $RLU_{dATP} \leq 10$ using a PhotonMaster or Lumitester C-110, you are below the low- detection limit.

NOTE: If $RLU_{dATP} \leq 50$ using a PhotonMaster or Lumitester C-110, consider accounting for background (RLU_{bg}). See Test Kit Instructions for guidance.

Dissolved ATP (dATP) Calculation:

$$dATP (pg \text{ ATP/mL}) = \frac{RLU_{dATP}}{RLU_{ATP1}} \times 10,000 (pg \text{ ATP / mL})$$

Calculations

NOTE: When using the QG21I tATP – only kit, skip final calculations and interpret tATP results as you would cATP results using the Interpretation Guidelines.

NOTE: If the results show for a given sample that $dATP (ng/mL) > tATP (ng/mL)$, report $cATP^{TM} = 0$ and $BSI^{TM} = 100\%$

Cellular ATP (cATP) Calculation:

$$cATP \left(\frac{ng \text{ ATP}}{mL} \right) = tATP \left(\frac{ng \text{ ATP}}{mL} \right) - dATP \left(\frac{ng \text{ ATP}}{mL} \right)$$

Microbial Equivalent (ME/mL):

$$cATP \left(\frac{ME}{mL} \right) = cATP (pg \text{ ATP/mL}) \times \frac{1 \text{ ME}}{0.001 \text{ pg ATP}}$$

NOTE: 1 ME (Microbial Equivalent) assumes 0.001 pg (1 fg) ATP per cell.

Biomass Stress Index (BSI) Calculation:

$$BSI (\%) = \frac{dATP (pg \text{ ATP/mL})}{tATP (pg \text{ ATP/mL})} \times 100\%$$

Interpretations Guidelines

Application	Good Control (pg cATP/mL)	Preventative Action (pg cATP/mL)	Corrective Action (pg cATP/mL)
Treated Process Water (Cooling, Bottom Water, Oilfield Non-Oxidizing Biocides or Non-Chemical Treatment)	< 100	100 to 1,000	> 1,000
Papermaking Product Quality (Newsprint, Fine papers)	< 1,000	1,000 to 10,000	> 10,000
Papermaking Odor Control (Paperboard, Recycle Water)	< 10,000	10,000 to 100,000	> 100,000

For BSI (when applicable), it can generally be interpreted that good control is achieved at levels of 75% or above. Preventive action should be taken at levels between 50% and 75%, and corrective action should be taken at levels below 50%.

NOTE: Interpretation Guidelines provided for general guidance. For best results, establish your own baseline and control levels.