



Pipet Filler, Safety Bulb

Proizvod #:	1465100
EUR Cijena:	Kontaktirajte nas
Otprema unutar 2 tjedna	

Unique cone adapter allows universal pipet filling. Simply squeeze rubber bulb to charge for filling, place pipet into cone adapter and insert into sample. Release the bulb to fill. Remove bulb and cover pipet end with finger for dispensing. Easy and safe. Autoclavable.

Specifikacije

Što je u pretincu:	Includes autoclavable elastic wristband.
Volumen/ veličina paketa:	1 each

Što je u pretincu

Includes autoclavable elastic wristband.