



## E.coli and total coliforms, pk/25

Proizvod #:	2495525
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Presence/Absence Broth with MUG ampules provide an easy-to-use test, ideal for screening drinking water samples for the presence of Total Coliforms and E. coli bacteria. Simply combine the P/A broth with your 100 mL sample and incubate for 24 to 48 hours at 35 °C. A colour change to yellow or yellow brown indicates presence of coliform bacteria.

**Complete Solutions for Testing Recreational Water**

**Complete Solutions for Testing Wastewater**

**Complete Solutions for Testing Drinking Water**

---

### Specifikacije

Materijal:	Glass
Volumen/ veličina paketa:	25 /pk