



## Ampules, m-ColiBlue24 broth, glass , pk/20

**Proizvod #:** 2608420  
**EUR Cijena:** Kontaktirajte nas  
**Nije prikazan datum isporuke**

**Important Note:** m-ColiBlue24 should only be used with EO sterilized pads that are currently manufactured by the following companies/ brands: Pall Corporation (e.g., Sterile Petri Dishes with Absorbant Pads), Private Label/OEM products manufactured by Pall Gelman, Sartorius (cellulose and glass fiber pads), Sartorius BioSart products with the integrated Petri dish, and Millipore EO sterilized pads.

No other pad is compatible with m-ColiBlue-24, and the use of any non-EO pad or EO pad, other than those manufactures recommended by Hach may result in unreliable and inaccurate testing.

**Read and confirm results in 24 hours**

**Get superior sensitivity - 1 CFU/100mL**

**Conduct emergency testing and routine monitoring**

**Eliminate special equipment and ultraviolet lamp**

**Minimize background growth of non-coliform bacteria**

---

### Specifikacije

Materijal:	Glass
Media Type:	m-ColiBlue24 Broth
Opis:	<p>In 24 hours, Hach's m-ColiBlue24® Broth give you simultaneous results for total coliforms and E. coli. Hach's ready-to-use media eliminates measuring, mixing, and autoclaving steps necessary to prepare dehydrated media. Glass ampules (26084-20), plastic ampules (260845-00), 100mL bottles (26084-42) of m-ColiBlue24 and agar plates (28052-15), afford maximum shelf life and ease of use for pre-measured and prepared media. Enumerate total coliforms and E. coli on one petri dish. Differentiate colonies easily - red and blue indicate total coliforms and blue specifies E. coli. Get optimal recovery of stressed and injured organisms.</p> <p>Use m-ColiBlue24 Broth to accurately monitor and evaluate total coliform and E. coli in drinking water; wastewater; bottled water; beverages; surface, ground, and well water; and ultrapure, chemical processing and pharmaceutical processing waters.</p>
Volumen/ veličina paketa:	20 /pk