



Aplikacijski softver za analizu u pivovarama

Proizvod #:	LZW942
EUR Cijena:	Kontaktirajte nas
Dostupno	

The application software LZW942 contains 12 brewery assays that conveniently upload to a DR6000 via USB. Assays are based on published and observed brewing methods according to MEBAC and include procedures and calculations for: Anthocyanogens, Iron, Steam volatile phenols, Bitter units, Photometric iodine, Thiobarbituric acid number (TAN), Free amino nitrogen, Reductones, Vicinal diketones.

Methods are ready to go.

The methods are convenient to upload via USB.

Keep control of time depending multiple measurements.

A timer function alerts about the next step and starts evaluation of multiple samples automatically.

Get the final result immediately after the analysis.

The photometer does the complex result calculation for you.