



Sodium Periodate Powder Pillows, 10 mL, pk/100

Proizvod #:

2107769

EUR Cijena:

Kontaktirajte nas

Dostupno

Opasno po zdravlje



Powder reagent in sealed packages for the determination of Manganese, cold periodate method. For 10 mL sample.

For determination of manganese by the periodate oxidation method. EPA approved for reporting wastewater analyses (digestion required). Hach method 8034. Sample size: 10 mL Pack of 100 Powder Pillows.

Specifikacije

Broj testova:	100
Naziv metode:	Periodate Oxidation
Platforma:	Praškasti reagensi PP
Volumen/ veličina paketa:	100 /pk