



Manganese Reagent Set, 0.006-0.700 mg/L Mn, 25 mL

Proizvod #:
EUR Cijena:
Otprema unutar 2 tjedna

2243300
Kontaktirajte nas

Opasno po zdravlje



For the determination of Manganese with PAN Method. For 25 mL sample.

For determination of Manganese by the PAN Method. Hach Method 8149. Range: 0.006 - 0.700 mg/L Mn. Sample Size: 25 mL. Reagent Set includes: Ascorbic Acid Powder Pillows, Alkaline-Cyanide Reagent and PAN Indicator Solution (0.1%) for 100 tests.

(DR/3000, DR/2000)

Specifikacije

Broj metode:	8149
Broj testova:	100
Mjerni rang:	0,006 - 0,700 mg/L Mn
Naziv metode:	PAN
Posebna značajka:	Manganese Reagent set, 25ml
Uređaj:	DR/3000, DR/2000
Volumen uzorka:	25 mL