



## Sodium periodate powder pillows for manganese, 25 mL, pk/100

**Proizvod #:**

**98499**

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Opasno po zdravlje



Powder reagent in sealed package for the determination of Manganese, cold periodate method, and Iron titrations. For 25 ml sample.

For determination of manganese by the periodate oxidation method. EPA approved for reporting wastewater analyses (digestion required). Hach method 8034. Sample size: 25 mL Pack of 100 Powder Pillows.

---

### Specifikacije

|                           |                       |
|---------------------------|-----------------------|
| Broj testova:             | 100                   |
| Naziv metode:             | Periodate Oxidation   |
| Platforma:                | Praškasti reagensi PP |
| Volumen/ veličina paketa: | 100 /pk               |