



Komplet reagensa za tanin i lignin

Proizvod #:

2244600

EUR Cijena:

Kontaktirajte nas

Dostupno

Opasno po zdravlje



For determination of Tannin and Lignin by the Tyrosine Method. For 25 mL sample.

Za određivanje tanina i lignina Tyrosine metodom. Hach metoda 8193. Raspon: 0,1-9,0 mg/L. Veličina uzorka: 25 mL. Komplet uključuje: otopinu natrijevog karbonata i TanniVer 3 reagens za tanin i lignin za do100 testova.

Specifikacije

Broj metode:	8193
Broj testova:	100
Mjerni rang:	0,1 - 9,0 mg/L kao taninska kiselina
Naziv metode:	Tyrosine
Parametar:	Tanin i lignin
Platforma:	Otopina
Uređaj:	Sve osim DR/820, PC II