



Reagent powder pillows, ChromaVer 3, chromium, 5 or 10 mL, pk/1000

Proizvod #:

1271028

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



For determination of hexavalent chromium by the Diphenylcarbohydrazide method. USEPA accepted for reporting for wastewater analyses using Hach Method 8023, Diphenylcarbohydrazide Method. Photometric Sample Size: 10 mL. Photometric Range: up to 0.70 mg/L. Kit Sample Size: 5 mL. Kit Range: Up to 1.50 Pack of 1000 powder pillows.

Specifikacije

Broj metode:	8023
Broj testova:	1000
EPA compliant:	Da
Mjerni rang:	0,010 - 0,700 mg/L Cr ⁶⁺
Naziv metode:	1,5-Diphenylcarbohydrazide
Parametar:	Chromium
Platforma:	Praškasti reagensi PP
Uređaj:	All except DR/3000, DR/2000, DR 900