



CDTA Digital titrator cartridge, 0.0800 M

Proizvod #: 1440201
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

Cyclohexanediamine tetraacetic acid standard solution. For Calcium and Total Hardness determination by Digital Titrator titration. 0.0800 M. Digits x 0.1 = mg/L CaCO_3 in 100 mL. Digits x 0.0001 = mmole.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikacije

Mjerni rang:	10 - 160 mg/L CaCO_3
Naziv metode:	CDTA Titration
Platforma:	Digitalni Titrator

Potrebni dodaci

- Zaštitne rukavice veličine L, plave, nitril, bez pudera, 50 komada (Item SM743L)