



Manganese Reagent Set, Powder Pillows, periodate oxidation, 0.1-20.0 mg/L Mn, 25 mL sample

Proizvod #:

2243200

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Opasno po zdravlje



For determination of Manganese by the Periodate Oxidation Method, Hach Method 8034. EPA Approved for wastewater reporting (digestion required) - Federal Register, 44, (116) 34193 (June 14, 1979). Range: 0.1 - 20.0 mg/L Mn. Sample Size: 25 mL. Set Contains: Buffer Powder Pillows (Citrate type) and Sodium Periodate Powder Pillows for 100 tests.

Specifikacije

Broj metode:	8034
Broj testova:	100
EPA compliant:	Da
Mjerni rang:	0,1 - 20,0 mg/L Mn
Naziv metode:	Periodate Oxidation
Parametar:	Manganese
Platforma:	Praškasti reagensi PP
Uređaj:	DR/3000, DR/2000
Volumen uzorka:	25 mL