



Be Right™



Amino acid F dilution solvent, for rapid liquid silica, 475 mL

Proizvod #:

2353011

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



For determination of Ultra Low Range silica by the Heteropoly Blue Rapid Liquid method. Hach method 8282 Range: 3 - 1000 µg/L SiO₂. Used with Amino Acid F Reagent Powder (2651155). Rapid Liquid™ methods offer an economical way to analyze multiple samples and to measure trace concentrations. Requires Pour-Thru Cell.

Specifikacije

Broj testova:	250
Naziv metode:	Heteropoly Blue
Platforma:	Otopina

Potrebni dodaci

- Zaštitne rukavice veličine L, plave, nitril, bez pudera, 50 komada (Item SM743L)