



Reagent set, phosphorus, ascorbic acid rapid liquid[™]

Proizvod #:

2678600

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Opasno po zdravlje



For determination of reactive Phosphorus (ortho-phosphate) by the Ascorbic Acid Rapid Liquid method. Hach Method 10055. Range: 19 - 3000 µg/L PO₄. Set includes: Ascorbic Acid Reagent Dilution Solution, Ascorbic Acid Reagent Powder and Molybdate Reagent Solution for approximately 250 tests.

Convenience

Rapid Liquid methods utilize liquid reagents, bottletop dispensers and pour-thru cells to provide simplified testing.

Performance

Since the pour-thru cell provides a stationary optical system, Rapid Liquid methods typically provide better precision than methods that require movement of sample cells.

Efficiency

Rapid Liquid methods are ideal for high-volume testing.

Specifikacije

Broj metode:	10055
Broj testova:	250
Mjerni rang:	19 - 3000 µg/L PO ₄
Naziv metode:	Ascorbic Acid
Platforma:	Otopina
Uređaj:	All except DR 900, DR/800 Series, PC II