



Free and total chlorine test strips, 0 - 10 mg/L, 50 tests

Proizvod #: 2745050
EUR Cijena: Kontaktirajte nas
Otprema unutar 2 tjedna

In the field, technicians can test many samples in only a few minutes, and make immediate evaluations on-site. No measuring, set-up, clean-up, or chemical handling are necessary. Lab technicians also use the test strips for pre-test screening-to detect the presence of materials that might interfere with testing. Steps: (mg/L) 0, 0.5, 1.0, 2.0, 4.0, 10.0.

Specifikacije

Jedinica:	50 tests
Mjerni rang:	0 - 10 mg/L Cl_2
Najmanji korak porasta:	Steps: 0, 0.5, 1, 2, 4, 10
Napomena:	*mg/L unless otherwise noted; ppb = $\mu\text{g/L}$; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Parametar:	Free and Total Chlorine
Platforma:	Testne trakice-Test Strip
Test Kit Model:	TS