



Filter, glass fiber, diameter 55 mm, preweighed standard methods TSS procedure, 100/pk

Proizvod #: 2546120
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

Specifikacije

Jedinica:	100/pk
Materijal:	Glass Microfiber
Promjer:	55 mm
Veličina pora:	1.5 μ m
Volumen/ veličina paketa:	100 /pk