



Sampling bag, sterile with dechlorinating agent, 100 mL, 100 pcs

Proizvod #: 2075333
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Contains non-nutritive pill with 10mg sodium thiosulfate (will dechlorinate up to 10 mg Cl_2 in 100 mL sample) four 100mL fill-lines.

Specifikacije

Kapacitet:	177 mL ml Fill Line
Volumen/ veličina paketa:	100