



Digital Titrator with plastic case

Proizvod #:	1690001
EUR Cijena:	Kontaktirajte nas
Otprema unutar 2 tjedna	

Portable and robust digital Titrator for fast and precise titrations. High resolution and high accuracy. Easy handling. Applications for more than 40 parameters with interchangeable cartridges.

Incl. plastic case, manual & accessories.

The Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field. Interchangeable Cartridges allow to change titrant in seconds. No cleaning or rinsing of burets is necessary. The number used for calculation of sample concentration is recorded directly from digital counter. For many methods, the displayed value represents the concentration directly.

Accurate results

Fast and Flexible

Versatile

Specifikacije

Dodaci uključeni:	Ne
Garancija:	24 mjeseci
Masa:	0,132 kg
Točnost:	± 1 % (Uncertainty of readings is 1 digit. Most samples require more than 100 digits)
Uključuje :	Digital Titrator only