



Test Kit, chlorine (free & total) & pH, model CN-67

Proizvod #:

1411100

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Opasno po zdravlje



Control of chlorine and pH in pools and spas to maintain safe bathing conditions and to protect against system corrosion. Supplied with color discs and comparator for monitoring chlorine and pH, the Free & Total Chlorine and pH Test Kit provides more accuracy than test kits using color cubes or step-wise comparators.

One of Hach's most popular test kits

Includes color discs and a comparator for monitoring chlorine and pH

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Broj testova:	300
Mjerni rang:	6.6 - 8.4 pH units
Mjerni rang (2):	0 - 3,4 mg/L
Model:	CN-67
Parametar:	Chlorine (Free & Total) and pH
Platforma:	Disk sa bojama-Color Disc
Primjena:	Pool and Spa