



Bottle, BOD, 300 mL, 24/pk

Proizvod #:	62124
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

The standard for BOD analysis. Borosilicate clear glass bottles with flared lip and ground-glass stoppers to prevent air entrapment and for water seal. Frosted marking area. Stoppers included. Pack of 24.

Specifikacije

Body Material:	Glass
Kapacitet:	300 mL
Materijal:	Borosilicate Glass
Package size UOM:	24pk
Volument/ veličina paketa:	24 /pk