



Beaker, PP, low form, 100 mL

Proizvod #:	108042
EUR Cijena:	Kontaktirajte nas
Dostupno	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Bull nose lip and molded pour spout minimise spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifikacije

Jedinica:	1 each
Kapacitet:	100 mL
Materijal:	Polypropylene, Low Form
Materijal tijela mjerača:	Polypropylene
Promjer, vanjski:	63.0 mm with spout
Volumen/ veličina paketa:	1 each