



Beaker, polypropylene, low form, 400 mL, 6/pk

Proizvod #:	108078
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifikacije

Jedinica:	6 /pk
Kapacitet:	400 mL
Materijal:	Polypropylene, Low Form
Materijal tijela mjerača:	Polypropylene
Promjer, vanjski:	99.1 mm with spout