



Be Right™



Flask, Filtering, Polypropylene, 1000 mL

Proizvod #:	2089953
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Polypropylene filtering flask with side-arm at 45° angle prevents tipping. Barbed side-arm accepts 1/4-3/8" I.D. tubing for application of vacuum. See page 380 for tubing; see page 375 for stopper.

Specifikacije

Jedinica:	1 each
Kapacitet:	1000 mL
Materijal tijela mjerača:	Polypropylene