



Water, deionised and demineralised, 500mL

Proizvod #:	27249
EUR Cijena:	Kontaktirajte nas
Dostupno	

Demineralised Type 2 Water

Specific Conductance: $\leq 1 \mu\text{S}/\text{cm}$,

Silica (SiO_2): $\leq 30 \mu\text{g}/\text{L}$,

Chloride (Cl^-): $\leq 20 \mu\text{g}/\text{L}$,

Hardness: Negative,

Consumption of KMnO_4 : Negative

Specifikacije

Opis:	Water, Deionized
Platforma:	Otopina
Volumen:	500 mL