



## Otopina za punjenje, referentna, 3 M KCl, 125 mL

|             |                   |
|-------------|-------------------|
| Proizvod #: | LZW9510.99        |
| EUR Cijena: | Kontaktirajte nas |
| Dostupno    |                   |

Otopina za punjenje referentnog dijela elektrode, kalijev klorid (KCl), 3 M, boca od 125 mL

---

### Specifikacije

|                           |                                               |
|---------------------------|-----------------------------------------------|
| Formula:                  | KCl 3M                                        |
| Opis:                     | PHC805, 5014T, 5021T, 5057, 5208, 5209, 5015T |
| Platforma:                | Otopina                                       |
| Solution Type:            | Ponovno punjenje                              |
| Volumen/ veličina paketa: | 125 mL                                        |