



KF1000 Series Volumetric Karl Fischer Titrator with 1 Burette and 2 Pumps

Proizvod #: KF1121.98

EUR Cijena: Kontaktirajte nas

Otprema unutar 3 dana

Simple and accurate volumetric titrations to measure moisture content

Simplify complex volumetric Karl Fischer titration with the Titrilab KF1000 from Hach. Straight from the box the KF1000 to easily and quickly perform moisture content determinations. The automatic system gives you more reliable results without complicated analysis. To make setup and titration simple, the KF1000 utilizes the application specific kit to make it quick and easy for anyone to set up and operate a test.

Specifikacije

Automatsko dodavanje reagensa:	2 integrated peristaltic pumps for reagent addition
Bireta: volumeni dostupni:	2.5/5/10/25 mL, ISO8655-3
Burette for Titrant:	For Titrant 1
Certificiranja:	Safety IEC/EN 61010-1; EMC IEC/EN 61326-1
Čuvanje uređaja:	-5 - 40 °C
Elektroda:	MTC Intellical ("Plug & Play" Digital)
Garancija:	24 mjeseci
Intervali/uzbune/podsjetnik kalibracije:	Da, korisnički definirano razdoblje
Izvoz podataka:	On USB flash drive, CSV format,
	Excel compatible
Kalibracija:	Titrant
Masa:	4 kg
Metrijal:	PP Latene, silicone, hardened glass, POM, stainless steel. Splashproof design.
Mjerna metoda:	Volumetric Karl Fischer, AC or DC imposed current
Mjerni načini:	Sample, blank, sample with blank, QC sample, QC sample with blank
Model:	KF1121
Pohrana podataka/memorija:	Last 100 samples, QC & blank analyses, last 10 calibrations
Preporučeno za:	AP0014.KF1121 Moisture Content
Prikaz krivulja kalibracije:	KF Titration curve
Primjena:	5
Pump flow rate:	170 mL/min
Pumpe:	2

Radna temperatura:	15 - 35 °C
Radna vlažnost:	20 - 80 %
Stalak za uzorak:	Integrated, magnetic stirring, KF glass beakers 150 mL
Što je u pretincu:	1 ea. KF1121-Titralab KF1000 Titrator
	1 ea. LZE144- Power Supply
	1 ea. LZE141- 4 position tube holder
	1 ea. Basic User Manual
	Sold Separately: Don't forget to order your Application Pack and Reagents and consider purchasing an AS1000 Sample Changer.
Supported peripherals:	Printer, PC software, balance, keyboard, mouse, barcode reader
Volumeni čaše prihvaćeni:	System specific beakers: Proprietary glass measurement cell, screw-on
Vrsta:	Karl Fischer
Zahtjevi napajanja (Hz):	50/60 Hz
Zahtjevi napajanja (voltaža):	100/240 V AC
Zaštita lozinkom:	Yes, user defined

Što je u pretincu

1 ea. KF1121-Titralab KF1000 Titrator 1 ea. LZE144- Power Supply 1 ea. LZE141- 4 position tube holder 1 ea. Basic User ManualSold Separately: Don't forget to order your Application Pack and Reagents and consider purchasing an AS1000 Sample Changer.

Potrebni dodaci

- Titration Application Pack Karl Fischer: Moisture Content (Item AP0014.KF1121)