



Kivetni test za sulfide 0,1-2,0 mg/L S²⁻

Proizvod #:

LCW053

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



For determination of Sulphide by the Dimethyl-p-phenylenediamine method.

Specifikacije

Broj testova:	35
Čuvanje uređaja:	15 °C - 25 °C
Mjerni rang:	0,1 - 2,0 mg/L S ²⁻
Naziv metode:	Dimetil-p-fenilendiamin
Opis:	Komplet reagensa za sulfide
Parametar:	Sulfid