



pH electrode for general use, 10 m fixed cable

Proizvod #: LZX472
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

Suitable for pH measurements in waste water treatment plants. Without in-built temperature element.

Specifikacije

Dijafagma:	Double ring diaphragm, porous PTFE
Dubina uranjanja:	105 mm
Duljina:	140 mm
Duljina kabla:	10 m fixed cable
Elektrolit:	KCl/AgCl + KNO ₃ Gel
Garancija:	6 mjeseci
Konektori:	Fixed cable
Mjerni rang:	0 - 14 pH
Primjena:	Wastewater
Radna temperatura:	0 - 80 °C
Rang tlaka:	3,5 bar
Reference system:	Ag/AgCl
Referentna elektroda:	Ag/AgCl
Temperaturni senzor:	None