



Beaker, Polypropylene, Low Form, 250 mL, 6/pk

Proizvod #:	108076
EUR Cijena:	Kontaktirajte nas
Otprema unutar 2 tjedna	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifikacije

Jedinica:	6 /pk
Kapacitet:	250 mL
Materijal:	Polypropylene, Low Form
Materijal tijela mjerača:	Polypropylene
Promjer, vanjski:	83.6 mm with spout