



**Be Right™**



## **Funnel, analytical, 254 mL approx. Volume**

**Proizvod #:**

**108370**

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Standard long stem funnel meets the needs of a wide variety of applications. Mounded ribbing prevents "air lock" commonly occurring during funneling and filtrations. A steep, 58-degree cone slope makes for fast and easy flow of either liquids or small solids. Polypropylene construction is resistant to chemicals and makes this funnel fully autoclavable for sterilization.

---

### **Specifikacije**

Volument/ veličina paketa:

1 each