



## Buret, Acrylic, 50 mL

<b>Proizvod #:</b>	<b>1405941</b>
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Durable, crystal-clear acrylic body is individually calibrated to meet the accuracy requirements of ASTM E287. Leak-proof stopcock with self-lubricating PTFE TFE plug. Excellent for general lab or field use in schools, industrial, municipal or food applications. Unaffected by dilute mineral acids and bases (except hydrofluoric acid and ammonium hydroxide). Can contain liquids within a pH range of 1 - 14 and up to 5% wt/wt of acids or 30% bases can be used. Not designed for exposure to alcohol or organic solvents.

---

### Specifikacije

Graduations:	0.1 mL
Kapacitet:	50 mL
Materijal tijela mjerača:	Acrylic
Posebna značajka:	Buret