



Tube, culture, 16 x 150 mm, screw cap, 10/pk

Proizvod #:	1436120
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Cap has a firmly bonded rubber liner providing a sure and highly inert sealing surface.

Specifikacije

Dimenzije (D x D):	16 mm x 150 mm
Duljina:	150 mm
Kapacitet:	Capacity to Neck 20 mL
Posebna značajka:	Tube, 16 x 150 mm
Promjer:	16 mm
Volumen/ veličina paketa:	10 /pk