



Iron and manganese color disc test kit, model IR-20

Proizvod #:

146300

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Opasno po zdravlje



When customers complain of stained and dingy laundry, this kit lets you do quick, on-the-spot tests for iron and manganese to see where the problem lies.

Method / Chemistry: Color disc colorimetry. Test: Iron/ Manganese. Range (mg/L): 0-5/ 0-3.

Simple, fast, color disc method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Mjerni rang:	0 - 3,0
Mjerni rang (2):	0 - 4,0 mg/L Fe
Model:	IR-20
Najmanji korak porasta:	0.2
Parametar:	Iron and Manganese
Platforma:	Disk sa bojama-Color Disc
Primjena:	Water Quality