



B4M TOC Analyser 6 Month Service Kit

Proizvod #:

19-KIT-125

EUR Cijena:

Kontaktirajte nas

Dostupno

Opasno po zdravlje



For BioTector B7000 MR analysers supplied from July 2014 to August 2015 with

analyser serial number format BAC-4xxx-xxxx

Specifikacije

Što je u pretincu:

Includes: Omron Relays (3), EMPP Tubing, Y Tube fitting (2), SR25 pump tube rail (2), PTFE Diaphragm, CO₂ Filter, Ozone Destructor catalyst, Filter Mats for vent (2), Lens Tissue for CO₂ Analyser, FKM/FPM o-ring (2)

Što je u pretincu

Includes: Omron Relays (3), EMPP Tubing, Y Tube fitting (2), SR25 pump tube rail (2), PTFE Diaphragm, CO₂ Filter, Ozone Destructor catalyst, Filter Mats for vent (2), Lens Tissue for CO₂ Analyser, FKM/FPM o-ring (2)