



Bottle, BOD, 60 mL capacity, marked 30 mL graduation line

Proizvod #:	190902
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included (60 mL sizes).

Specifikacije

Body Material:	Borosilicate Glass
Kapacitet:	60 mL
Materijal:	Borosilicate Glass
Package size UOM:	1pk
Volumen/ veličina paketa:	1 each