



Tubing, nalgene 8000, vacuum, 9.5 mm (I.D.)

Proizvod #:	2074146
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Similar to Nalgene 8000 regular tubing, but with heavier walls. Imprinted in red for instant identification. Saves time and avoids errors which could invalidate experiments or injure laboratory personnel.

Specifikacije

Debljina stijenke cijevi:	1/4 in
Duljina:	50 ft
Formula (Text):	R-3603
I. D.:	3/8 in
Materijal tijela mjerača:	Tygon
Promjer, vanjski:	7/8 in