



Bag, sterile, WhirlPak w/ dechlorinating agent, 100 mL, 25/ pk

Proizvod #:	2075325
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Specially designed for chlorinated water samples (Complies with USEPA and Standard Methods criteria). Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Contains non-nutritive pill with 10 mg sodium thiosulfate (will dechlorinate up to 10 mg/ Cl_2 in 100 mL sample) Four ounce and 100 mL fill-lines. "Write-on" strip accepts waterproof ink.

Specifikacije

Kapacitet:	177 mL
Volumen/ veličina paketa:	25