



Cylinder, graduated, mixing, Glass, 50 mL ± 0.4 mL, 1.0 mL divisions, glass stopper #16

Proizvod #: 2088641
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

This tall form glass graduated cylinder is calibrated “to contain” with a molded in base. Each includes glass penny head stopper. Graduations are in durable white enamel bake.

Specifikacije

Certificiranja:	Class B
Kapacitet:	50 mL
Materijal tijela mjerača:	Glass
Mjerna metoda:	To Contain
Volumen/ veličina paketa:	1 each
Vrsta:	Jednometrički, miješanje