



Pipet rack

Proizvod #:

2094700

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Polypropylene, fifty-place rack for safe, orderly storage of pipets. Helps create an efficient storage system or a well organized liquid handling station. Holes in top and middle decks are 16 mm diameter (0.63"); bottom holes are tapered to 6 mm (0.24") to prevent contact with supporting surface. Autoclavable to 122°C. Overall dimension: 21 x 11 x 22 cm (8.3 x 4.3 x 8.7").

Specifikacije

Dubina:	8.7 Inches
Duljina:	210 mm
Širina:	4.3 Inches
Volumen/ veličina paketa:	1 each