



Crucible, Gooch, 40 mL Capacity

Proizvod #:	2098641
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Perforated bottom, glazed inside and out except for outside bottom surface. (Perforations are 1 mm in diameter). Filter for Gooch Crucible available and sold separately.

25 mL Gooch Crucible (# 2098640) is recommended for use in Standard Methods 2540 C and 2540 D to determine Total Suspended Solids.

Specifikacije

Kapacitet:	40 mL
Volumen/ veličina paketa:	1 /pk