



Bag, sterile, WhirlPak, 207 mL, 100/pk

Proizvod #:	2233199
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Sampling and Handling

Whirl-Pak Sterile Bag Vol: 207 mL (7oz), Dim (W x L): 9.5 cm x 18 cm, (3.75 x 7") 3.0 mil thick

Complies with Standard Methods and USEPA for drinking water sample collection. Made from poly resins approved by USDA and FDA for food contact. Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Suitable for liquid or solid sample collection. "Write-on" surface accepts waterproof ink. 24 oz bag perfect for incubating petri dishes in water bath applications.

Specifikacije

Kapacitet:	207 mL
Volumen/ veličina paketa:	100