



Bottle, Sampling, Polycarbonate, 1000 mL, 4/pk

Proizvod #:	2324663
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Transparent polycarbonate bottle with cap features wide mouth (3.8 cm/1.48" diameter) opening for easy filling and square design for maximum space efficiency. Autoclavable with molded-in graduations (do not engage threads for autoclave cycle.)

Specifikacije

Body Material:	Polycarbonate
Kapacitet:	1000 mL
Volumen/ veličina paketa:	4 /pk