



Filter holder, for vacuum filtration, glass

Proizvod #:	234000
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Beaker-style funnel (300 mL) with 47 mm membrane support for filtration of aqueous samples. Fritted borosilicate glass membrane support for even filtration across surface of membrane. Funnel attaches evenly to base via aluminum clamp. Filter flask sold separately.

Specifikacije

Promjer:	47 mm
Što je u pretincu:	Borosilicate glass funnel and base; anodized aluminum clamp; support screen (sintered glass); rubber stopper (No. 8).
Volumen/ veličina paketa:	1

Što je u pretincu

Borosilicate glass funnel and base; anodized aluminum clamp; support screen (sintered glass); rubber stopper (No. 8).