



Rack, Test Tube, Polyethylene, 13 mm (O.D.), 90 Holes

Proizvod #:	2497900
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Ninety-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

Specifikacije

Aranžman testova:	6 x 15
Promjer:	13 mm
Rupe za testove:	90
Volumen/ veličina paketa:	1 each