



Test tube/ MPN rack, polyethylene, 20 mm (O.D.), 40 holes

Proizvod #:	2497902
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Forty-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

Specifikacije

Aranžman testova:	4 x 10
Promjer:	20 mm
Rupe za testove:	40
Volumen/ veličina paketa:	1 each