



Bottle rack, plastic-coated, 20 compartments

Proizvod #:	2551801
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Plastic-coated wire rack provides corrosion resistant transport system for sample bottles and bags. Racks feature 15 or 20 compartments with double-grid bottom to prevent loss of samples during transport. 15 compartment rack is ideal for Whirl-Pak sample bags (4 oz). 20 compartment rack fits most bottles associated with standard drinking water sampling. Autoclavable.

Specifikacije

Posebna značajka:	Rack, 20 Bottle Compartments
Volumen/ veličina paketa:	1