



pH Paper, 0 - 14 pH Range, 100 tests

Proizvod #:	2601300
EUR Cijena:	Kontaktirajte nas
Dostupno	

Provides a simple, fast, economical method of pH measurement. Obtain accurate results in weakly buffered, turbid and lightly colored samples. Measures 0 - 14 pH.

Test Strips are great for quick evaluations of pH levels in the lab or in the field

Easy to use, disposable, and inexpensive!

Specifikacije

Broj metode:	10056
Broj testova:	100
Jedinica:	100 tests
Mjerni rang:	0 - 14 pH units
Najmanji korak porasta:	0, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14
Parametar:	pH
Platforma:	Testne trakice-Test Strip