



Bottle, wash, diamond alert, 500 mL, Variety 6/pk

Proizvod #:	2613049
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Wide-mouthed and color-coded wash bottles provide safe and positive identification of common laboratory solutions. Standard Diamond Alert symbols indicate type of hazard contents may represent. Polyethylene bottles with color-coded polypropylene screw caps and spray-tip assembly. Variety 6-pack includes one bottle each for water, acetone, methanol, isopropanol, toluene, and ethanol.

Specifikacije

Boja:	Transparent
Jedinica:	6 pk
Kapacitet:	500 mL
Materijal:	LDPE
Stem:	Blue Cap/Stem
Volumen/ veličina paketa:	6 pk