



Sample vial, 120 mL w/ thiosulfate & custody tie, 200/pk

Proizvod #:	2634001
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Sample Vial with Thiosulfate and Custody Tie

120 mL sterile, hinged-cap vials meet all EPA requirements for microbiology sampling. Collection, incubation, and accurate reading can all be done in the vial without a transfer step. 100 mL raised fill line. Polypropylene. (**Please Note:** These sample vials are not autoclaveable)

Specifikacije

Duljina:	93mm
Jedinica:	200 /pk
Kapacitet:	fill line at 100 mL (+/-2.5 mL)
Posebna značajka:	Sample Vial, with Thiosulfate
Promjer:	43mm
Volumen/ veličina paketa:	200 /pk