



Buret, Class A, Glass, 50 mL

Proizvod #: 2636541
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

Glass, funnel-fill buret. Certified Class A (ASTM Specification E287) with permanent serial number and certificate for NIST traceability. PTFE stopcock does not require grease. Dust cap included. 10 mL buret meets requirements for potentiometric titrations (ASTM D664).

Specifikacije

Kapacitet:	50 mL
Materijal tijela mjerača:	Class A Glass
Posebna značajka:	Buret